

Product Datasheet

CDX2 antibody (orb387929)

Description

Mouse monoclonal antibody to CDX2


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunohistochemical
staining of human
Co...

Immunogen


Recombinant fragment (around aa150-249) of human CDX2 protein (exact sequence is proprietary)

Target

CDX2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.


Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Immunohistochemical
staining of human
Co...

Note

For research use only

Application notes

CDX2 (GI Epithelial Marker)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

40kDa

Uniprot ID

[Q99626](#)

Entrez

[1045](#)

Dilution Range

ELISA: 2-4 µg/ml, FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.