



Product Datasheet

CD66 antibody (orb387933)

Description	Mouse monoclonal antibody to CD66
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Human colon carcinoma extract
Target	CEACAM5
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Carcinoembryonic Antigen (CEA) / CD66
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	80-200kDa
Uniprot ID	P06731
Entrez	1048
Dilution Range	FACS: 0.5-1 ug/million cells, IF/ICC: 1-2 ug/ml, WB: 0.5-1 ug/ml, IHC-P: 0.5-1 ug/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Co...



Western blot analysis of human Stomach L...