

Product Datasheet

TNFSF13B antibody (orb387950)

Description	Mouse monoclonal antibody to TNFSF13B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant human CD257 protein
Target	TNFSF13B
Preservatives	200ug/ml of Ab Purified from rabbit anti-serum by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CD257 / BAFF / TNFSF13B
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	31kDa (Monomer)
Uniprot ID	Q9Y275
Entrez	10673
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml
Expiration Date	12 months from date of receipt.