

Product Datasheet

CFTR antibody (orb387952)

Description

Mouse monoclonal antibody to CFTR

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant human CFTR fragment

Target

CFTR

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

CFTR (Cystic Fibrosis Transmembrane Conductance Regulator)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

165-170kDa

Uniprot ID

P13569

Entrez


1080

Dilution Range


IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of mouse Ki...



Immunohistochemical staining of human Pa...