


Product Datasheet

CFTR antibody (orb387957)

Description	Rabbit monoclonal antibody to CFTR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant fragment (around aa 258-385) of human CFTR protein (exact sequence is proprietary)
Target	CFTR
Preservatives	200ug/ml of Ab Purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CFTR (Cystic Fibrosis Transmembrane Conductance Regulator)
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	165-170kDa
Uniprot ID	P13569
Entrez	1080
Dilution Range	IF/ICC: 0.5-1 µg/ml, IHC-P1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Pa...