

Product Datasheet

CGB antibody (orb387965)

Description	Mouse monoclonal antibody to CGB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Immunogen	Purified human HCG-beta
Target	CGB
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	HCG-beta (Pregnancy & Choriocarcinoma Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	22kDa
Uniprot ID	P01233
Entrez	1082
Dilution Range	IHC-Fr: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.