

Product Datasheet

CKB antibody (orb387986)

Description

Mouse monoclonal antibody to CKB

Species/Host

Mouse

Reactivity

Bovine, Canine, Frog, Gallus, Human, Monkey, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC-Fr

Applications
Immunogen

Human CKBB protein

Target

CKB

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Creatine Phosphokinase-BB (CK-BB)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

43kDa (Monomer); 86kDa (Dimer)

Uniprot ID

P12277

Entrez

1152

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-Fr: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.