

## Product Datasheet

### CCR5 antibody (orb387991)

**Description**

Mouse monoclonal antibody to CCR5

|                            |                                                                                                                                                                                  |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Species/Host</b>        | Mouse                                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                     |
| <b>Tested Applications</b> | FC, IF, IHC-Fr                                                                                                                                                                   |
| <b>Immunogen</b>           | Human native CCR5 protein                                                                                                                                                        |
| <b>Target</b>              | CCR5                                                                                                                                                                             |
| <b>Preservatives</b>       | 200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml. |
| <b>Concentration</b>       | Purified Ab with BSA and Azide at 200ug/ml                                                                                                                                       |
| <b>Storage</b>             | Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.                      |
| <b>Note</b>                | For research use only                                                                                                                                                            |
| <b>Application notes</b>   | CD195 (CC-Chemokine Receptor 5)                                                                                                                                                  |
| <b>Isotype</b>             | Mouse / IgG2a, kappa                                                                                                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                       |
| <b>MW</b>                  | 45kDa                                                                                                                                                                            |
| <b>Uniprot ID</b>          | <a href="#">P51681</a>                                                                                                                                                           |
| <b>Entrez</b>              | <a href="#">1234</a>                                                                                                                                                             |
| <b>Dilution Range</b>      | FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-Fr: 0.5-1.0 µg/ml                                                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                  |