

Product Datasheet

SERPINA3 antibody (orb388003)

Description

Mouse monoclonal antibody to SERPINA3


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of
SERPINA3...

Immunogen

Recombinant human Antichymotrypsin (AACT) protein fragment (around aa 49-187) (exact sequence is proprietary)

Target

SERPINA3

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Alpha-1-Antichymotrypsin (SERPINA3) (Histiocytoma Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

65-76kDa

Uniprot ID

[P01011](#)

Entrez

[12](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.