

Product Datasheet

CTNNB1 antibody (orb388013)

Description

Mouse monoclonal antibody to CTNNB1

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IF, WB

Applications
Immunogen

Recombinant fusion protein consisting of the N-terminal half of human beta-Catenin fused to maltose binding protein

Target

CTNNB1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -80°C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Catenin, beta (p120)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

92kDa

Uniprot ID
[P35222](#)
Entrez
[1499](#)
Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml,
WB: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.