

Product Datasheet

ENO2 antibody (orb388068)

Description

Mouse monoclonal antibody to ENO2

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

A synthetic peptide of human NSE gamma (around aa416-433) (exact sequence is proprietary)

Target

ENO2

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

NSE gamma (Neuron Specific Enolase, gamma) (Neuroendocrine Marker)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

~50kDa

Uniprot ID

P09104

Entrez


2026

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.1-0.2 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Ph...



Immunohistochemical staining of human Ce...



Immunohistochemical staining of mouse Pa...