

Product Datasheet

F13A1 antibody (orb388088)

Description

Mouse monoclonal antibody to F13A1


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC, WB

Applications

Immunohistochemical
staining of human
Pl...

Immunogen

Recombinant fragment of human Factor XIIIa protein (aa46-181) (exact sequence is proprietary)

Target

F13A1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Factor XIIIa (Coagulation Factor XIIIa Chain)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

83kDa

Uniprot ID

P00488

Entrez

2162

Dilution Range

ELISA: 2-4 µg/ml, FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.