


Product Datasheet

AMACR antibody (orb388129)

Description	Mouse monoclonal antibody to AMACR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human AMACR protein
Target	AMACR
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	AMACR / p504S (Prostate Cancer Marker)
Isotype	Mouse / IgG2b, lambda
Clonality	Monoclonal
MW	42kDa
Uniprot ID	Q9UHK6
Entrez	23600
Dilution Range	IF/ICC: 1-2 µg/ml, IHC-P1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Prostate Cancer tissue.