

Product Datasheet

TRIM29 antibody (orb388131)

Description

Mouse monoclonal antibody to TRIM29


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications

Immunohistochemical
staining of human
Lu...

Immunogen


Recombinant fragment (126 Amino acid residues
between aa 1-200) of human TRIM29 protein

Target

TRIM29

Preservatives

200ug/ml of Ab Purified from Bioreactor
Concentrate by Protein A/G. Prepared in 10mM PBS
with 0.05% rAlbumin & 0.05% azide. Also available
WITHOUT rAlbumin & azide at 1.0mg/ml.


Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Immunohistochemical
staining of human
Es...

Note

For research use only

Application notes

TRIM29 (Lung Squamous Cell Carcinoma Marker)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

66kDa

Uniprot ID

[Q14134](#)

Entrez

[23650](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-
P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.