



Product Datasheet

ALPL antibody (orb388137)


Description	Mouse monoclonal antibody to ALPL
Species/Host	Mouse
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human ALPL protein
Target	ALPL
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Alkaline Phosphatase (Tissue-Nonspecific)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	55kDa
Uniprot ID	P05186
Entrez	249
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human En...



Immunohistochemical staining of human Te...



Immunohistochemical staining of human Ov...