

Product Datasheet

FTL antibody (orb388144)

Description

Mouse monoclonal antibody to FTL

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

Recombinant fragment (around aa 38-165) of human FTL protein (exact sequence is proprietary)

Target

FTL

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Ferritin, Light Chain (FTL) (Microglia Marker)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

19-25kDa

Uniprot ID

[P02792](#)

Entrez


[2512](#)

Dilution Range


FACS: 0.1-0.2 µg/million cells, IF/ICC: 0.1-0.2 µg/ml, WB: 0.1-0.2 µg/ml, IHC-P: 0.1-0.2 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Pl...



Immunohistochemical staining of human Pa...



Western blot analysis of A431, HeLa, Liv...