


Product Datasheet

FUT4 antibody (orb388152)

Description	Mouse monoclonal antibody to FUT4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Myelomonocytic leukemia cells
Target	FUT4
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CD15 / FUT4 (Reed-Sternberg Cell Marker)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	~220kDa (Glycoprotein)
Uniprot ID	P51993, P22083
Entrez	2528, 10690, 2526, 2527, 2529
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ho...