

Product Datasheet

XRCC6 (phospho-70) antibody (orb388156)

Description Mouse monoclonal antibody to XRCC6 (phospho-70)

Species/Host Mouse

Reactivity Gallus, Human, Primate, Rat

Conjugation Unconjugated

Tested Applications FC, ICC, IF

Immunogen Nuclear extract of human HL-60 cells

Target XRCC6 (p70) & XRCC5 (p80)

Preservatives 200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration Purified Ab with BSA and Azide at 200ug/ml

Storage Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note For research use only

Application notes Ku (p70/p80) (Nuclear Marker)

Isotype Mouse / IgG1, kappa

Clonality Monoclonal


MW 70kDa & 80kDa

Uniprot ID **P13010**, **P12956**


Entrez **80**, **70**, **7520**, **2547**

Dilution Range FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.25-0.5 µg/ml

Expiration Date 12 months from date of receipt.



Flow
cytometric
analysis of
K562 Cells
u...



Flow
cytometric
analysis of
293T Cells
u...