

Product Datasheet

FBXO5 antibody (orb388157)

Description

Mouse monoclonal antibody to FBXO5

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

Recombinant fragment (203 amino acid residues between aa 1-250) of human EMI1 protein

Target

FBXO5

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

EMI1 (Early Mitotic Inhibitor-1)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

56kDa

Uniprot ID

[Q9UKT4](#)

Entrez


[26271](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 1-2 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Ov...



Western blot analysis of Hela, HepG2, HE...



Immunohistochemical staining of human Ov...