



Product Datasheet

GJB1 antibody (orb388167)

Description	Mouse monoclonal antibody to GJB1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human GJB1 protein
Target	GJB1
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Connexin 32 (Gap Junction Protein)
Isotype	Mouse / IgG
Clonality	Monoclonal
MW	27-32kDa
Uniprot ID	P08034
Entrez	2705
Dilution Range	IF/ICC: 1-2 µg/ml, FACS: 0.5-1 µg/million cells, WB: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Liver



Immunohistochemical staining of human Testis