


Product Datasheet

CD57 antibody (orb388172)

Description	Mouse monoclonal antibody to CD57
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Membrane antigen from HSB-2 cells
Target	B3GAT1
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	CD57 / B3GAT1 (Natural Killer Cell Marker)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	~110kDa
Uniprot ID	Q9P2W7
Entrez	27087
Dilution Range	IF/ICC: 0.5-1 µg/ml, IHC-P: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Sp...