

Product Datasheet

ABO antibody (orb388182)

Description

Mouse monoclonal antibody to ABO

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Mixture of erythrocytes of blood group A and glycoprotein fraction isolated from the saliva of secretors with blood group A
Target	ABO
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Blood Group Antigen A (CD173)
Isotype	Mouse / IgM, kappa
Clonality	Monoclonal
MW	Multiple
Uniprot ID	P16442
Entrez	28
Dilution Range	IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.