


Product Datasheet

ABO antibody (orb388187)

Description	Mouse monoclonal antibody to ABO
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Mucin isolated from an ovarian cyst fluid
Target	ABO
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Blood Group Antigen A (CD173)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	Multiple
Uniprot ID	P16442
Entrez	28
Dilution Range	IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Co...