

Product Datasheet

Galectin-13 antibody (orb388190)

Description

Mouse monoclonal antibody to Galectin-13


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IF, IHC-P

Applications

Immunohistochemical
staining of human
Pl...

Immunogen

Recombinant human Galectin-13 protein fragment
(aa23-134)

Target

LGALS13

Preservatives

200ug/ml of Ab Purified from Bioreactor
Concentrate by Protein A/G. Prepared in 10mM
PBS with 0.05% rAlbumin & 0.05% azide. Also
available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Note

For research use only

Application notes

Galectin-13 (GAL13) / Placental Protein 13 (PP13)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

16kDa (monomer); 32kDa (homodimer)

Uniprot ID
[Q9UHV8](#)
Entrez
[29124](#)
Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml,
IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.