


Product Datasheet

GYPA antibody (orb388199)

Description	Rabbit monoclonal antibody to GYPA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human GYPA protein
Target	GYPA
Preservatives	200ug/ml of Ab Purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Glycophorin A / CD235a (Erythrocyte Marker)
Isotype	Rabbit / IgG
Clonality	Monoclonal
MW	39kDa
Uniprot ID	P02724
Entrez	2993, 2994
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human an...