

Product Datasheet

HLA-A antibody (orb388210)

Description	Mouse monoclonal antibody to HLA-A
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Human peripheral blood leukocytes
Target	HLA-A; HLA-B; HLA-C
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	HLA-ABC (MHC I)
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	~41kDa
Uniprot ID	P30443, P30499, P01889
Entrez	3106, 3107, 3105
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.