

Product Datasheet

HLA-DRB1 antibody (orb388232)

Description

Mouse monoclonal antibody to HLA-DRB1

Species/Host

Mouse

Reactivity

Canine, Human, Monkey, Primate, Rodent

Conjugation

Unconjugated

Tested

FC, IF

Applications
Immunogen

Activated human peripheral blood mononuclear cells

Target

HLA-DRB1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C.
Antibody without azide - store at -20 to -80°C.
Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

HLA-DRB (MHC II)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

~28kDa (beta chain)

Uniprot ID

P01911

Entrez

3123

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.