

Product Datasheet

HLA-DRB1 antibody (orb388233)

Description

Mouse monoclonal antibody to HLA-DRB1

Species/Host

Mouse

Reactivity

Human, Monkey

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

Activated human peripheral blood mononuclear cells (LN-3 & HLA-DRB/1067)

Target

HLA-DRB1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

HLA-DRB (MHC II)

Isotype

Mouse / IgG2b, kappa

Clonality

Monoclonal

MW

28kDa (beta chain)

Uniprot ID

[P01911](#)

Entrez


[3123](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.25-0.5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human To...



Immunohistochemical staining of human Hi...