

Product Datasheet

FOXA1 antibody (orb388235)

Description

Mouse monoclonal antibody to FOXA1

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

Recombinant human FOXA1 protein fragment (aa372-472) (exact sequence is proprietary)

Target

FOXA1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

FOXA1 / HNF3A

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

49kDa

Uniprot ID

[P55317](#)

Entrez


[3169](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human prostate tissue



Immunohistochemical staining of human prostate tissue



Western blot analysis of MCF-7 (Lane 1) ...