


Product Datasheet

HSPG2 antibody (orb388249)

Description	Rat monoclonal antibody to HSPG2
Species/Host	Rat
Reactivity	Bovine, Fish, Human, Monkey, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Murine EHS laminin preparation
Target	HSPG2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Heparan Sulfate Proteoglycan (Large) / Perlecan
Isotype	Rat / IgG2a, kappa
Clonality	Monoclonal
MW	> 400kDa
Uniprot ID	P98160
Entrez	3339
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Te...