

Product Datasheet

TNC antibody (orb388251)

Description

Mouse monoclonal antibody to TNC

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of human
Lu...

Immunogen

Human breast carcinoma

Target

TNC

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate
by Protein A/G. Prepared in 10mM PBS with 0.05%
rAlbumin & 0.05% azide. Also available WITHOUT
rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Note

For research use only

Application notes

Tenascin C (Stromal Marker For Epithelial
Malignancy)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

210kDa and 300kDa

Uniprot ID

[P24821](#)

Entrez

[3371](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-
P: 2-4 µg/ml

Expiration Date

12 months from date of receipt.