

Product Datasheet

IGHA1 antibody (orb388265)

Description

Mouse monoclonal antibody to IGHA1

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Immunogen	Purified human IgA
Target	IGHA1 & IGHA2
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	IgA (Immunoglobulin Alpha Heavy Chain) (B-Cell Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	50-75kDa
Uniprot ID	P01877, P01876
Entrez	3493, 3494
Dilution Range	FACS: 1-2 µg/million cells, IF/ICC: 0.5-1.0 µg/ml, IHC-P: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.