



Product Datasheet

IGL antibody (orb388300)

Description	Mouse monoclonal antibody to IGL
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Purified human IgG (LcN-2 and ICO-106)
Target	IGL
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Lambda Light Chain (B-Cell Marker)
Isotype	Mouse / IgG1, kappa + Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	~22.5kDa
Uniprot ID	P01842, P01701
Entrez	3535, 3546
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1.0 µg/ml, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Intestina...



Immunohistochemical staining of human To...