



Product Datasheet


CYP3A43 antibody (orb38834)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	CYP3A43
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57670
Uniprot ID	Q9HB55
NCBI	NP_073731.1 , NP_476436.1 , NP_001265850.1 , NP_476437.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of 293 cell

...



Western blot analysis of in mouse brain ...