

Product Datasheet

KDR antibody (orb388344)

Description

Mouse monoclonal antibody to KDR

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies, IF
Immunogen	Recombinant human VEGF-R2 protein
Target	KDR
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	VEGFR2 / CD309 / Flk-1 / KDR3 (Hemangioblast Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	150-230kDa
Uniprot ID	P35968
Entrez	3791
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.