

Product Datasheet

ARG1 antibody (orb388346)

Description

Mouse monoclonal antibody to ARG1

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunohistochemical
staining of human
He...

Immunogen

Recombinant fragment (around aa11-97) of human ARG1 protein(exact sequence is proprietary)

Target

ARG1

Preservatives

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Arginase 1 (Hepatocellular Carcinoma Marker)

Isotype

Mouse / IgG3, kappa

Clonality

Monoclonal

MW

35-38kDa

Uniprot ID

P05089

Entrez

383

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 2-4 µg/ml

Expiration Date

12 months from date of receipt.