

Product Datasheet

KRT6A antibody (orb388349)

Description

Mouse monoclonal antibody to KRT6A


Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunohistochemical
staining of human
Bl...

Immunogen

A synthetic peptide of 11 amino acids (GSSTIKYTTTS) from C-terminus of human keratin 6

Target

KRT6A; KRT6B; KRT6C

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin 6 (KRT6) (Hyperproliferation-Related Keratin)

Isotype

Mouse / IgG2a, kappa

Clonality

Monoclonal

MW

56kDa

Uniprot ID

P04259, P02538, P48668

Entrez

3854, 3853, 286887

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.1-0.2 µg/ml

Expiration Date

12 months from date of receipt.