

Product Datasheet

KRT7 antibody (orb388357)

Description Mouse monoclonal antibody to KRT7

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, IHC-P

Immunogen Recombinant full-length human KRT7 protein

Target KRT7

Preservatives 200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration Purified Ab with BSA and Azide at 200ug/ml

Storage Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note For research use only

Application notes Cytokeratin 7 (KRT7) (Glandular and Transitional Epithelial Marker)

Isotype Mouse / IgG1, kappa + Mouse / IgG1, kappa


Clonality Monoclonal

MW 55kDa

Uniprot ID [P08729](#)
Entrez [3855](#)
Dilution Range FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date 12 months from date of receipt.


Immunohistochemical staining of human Ov...



Immunohistochemical staining of human Lu...



Immunohistochemical staining of human En...