

Product Datasheet

KRT7 antibody (orb388358)

Description

Mouse monoclonal antibody to KRT7

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Recombinant full-length human KRT7 protein (KRT7/760); OTN 11, ovarian carcinoma cell line (OV-TL12/30)

Target

KRT7

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin 7 (KRT7) (Glandular and Transitional Epithelial Marker)

Isotype

Mouse / IgG1, kappa + Mouse / IgG1, kappa

Clonality

Monoclonal

MW

55kDa

Uniprot ID

[P08729](#)

Entrez


[3855](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Ov...



Immunohistochemical staining of human Lu...