


Product Datasheet

KRT15 antibody (orb388385)

Description	Mouse monoclonal antibody to KRT15
Species/Host	Mouse
Reactivity	Bovine, Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human KRT15 protein
Target	KRT15
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Cytokeratin 15 (KRT15) (Basal Cell Marker)
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	52kDa
Uniprot ID	P19012
Entrez	3866
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-Fr: 0.5-1.0 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ba...