


Product Datasheet

KRT18 antibody (orb388391)

Description	Mouse monoclonal antibody to KRT18
Species/Host	Mouse
Reactivity	Bovine, Canine, Hamster, Human, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	PMC-42 human breast carcinoma cells
Target	KRT18
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Cytokeratin 18 (KRT18)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	45kDa
Uniprot ID	P05783
Entrez	3875
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela (Lane 1) A...



Immunohistochemical staining of human Co...