

Product Datasheet

KRT19 antibody (orb388409)

Description

Mouse monoclonal antibody to KRT19

Species/Host

Mouse

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

Recombinant full-length human KRT19 protein

Target

KRT19

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Cytokeratin 19 (KRT19) (Pancreatic Stem Cell Marker)

Isotype

Mouse / IgG2a, lambda + Mouse / IgG1, kappa

Clonality

Monoclonal

MW

40kDa

Uniprot ID

[P08727](#)

Entrez


[3880](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Ce...



Immunohistochemical staining of human Bl...



Immunohistochemical staining of human Ov...