

## Product Datasheet

### MCM7 antibody (orb388421)

**Description**

Mouse monoclonal antibody to MCM7

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

IHC-P, WB

**Applications**
**Immunogen**

Recombinant human MCM7 protein fragment (aa195-319)

**Target**

MCM7

**Preservatives**

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin &amp; 0.05% azide. Also available WITHOUT rAlbumin &amp; azide at 1.0mg/ml.

**Concentration**

Purified Ab with BSA and Azide at 200ug/ml

**Storage**

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

**Note**

For research use only

**Application notes**

MCM7 (Proliferation Marker)

**Isotype**

Mouse / IgG2b

**Clonality**

Monoclonal

**MW**

88kDa

**Uniprot ID**

[P33993](#)

**Entrez**


[4176](#)

**Dilution Range**

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

**Expiration Date**

12 months from date of receipt.



Immunohistochemical staining of human To...