

Product Datasheet

RPS6KB2 antibody (orb38844)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RPS6KB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53455
Uniprot ID	Q9UBS0
NCBI	NP_003943.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

