

Product Datasheet

ABCC1 antibody (orb388443)

Description

Mouse monoclonal antibody to ABCC1

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, WB

Applications

Immunohistochemical
staining of human
To...

Immunogen

Recombinant full-length human ABCC1 protein

Target

ABCC1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Immunohistochemical
staining of human
Pa...

Note

For research use only

Application notes

MRP1 / ABCC1 (Multidrug Resistance Related Protein 1)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

190kDa

Uniprot ID

[P33527](#)

Entrez

[4363](#)

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.5-1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.

