

Product Datasheet

MYC antibody (orb388485)

Description

Mouse monoclonal antibody to MYC

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

A synthetic peptide, corresponding to aa 408-439 (AEEQKLISEEDLLRKRREQLKHKLEQLRNSCA) from C-terminus of human c-myc, coupled to KLH.

Target

MYC

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

c-Myc Oncoprotein

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

62-64kDa

Uniprot ID

P01106

Entrez


4609

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Ce...