

Product Datasheet

MYH11 antibody (orb388491)

Description

Mouse monoclonal antibody to MYH11

Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Recombinant full-length human MYH11 protein

Target

MYH11

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Smooth Muscle Myosin Heavy Chain (SM-MHC) (Leiomyosarcoma & Myoepithelial Cell Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

205kDa (MHC-1) and 200kDa (MHC-2)

Uniprot ID

[P35749](#)

Entrez

[4629](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of human Br...



Immunohistochemical staining of human Le...



Immunohistochemical staining of human Co...