


Product Datasheet

MYH11 antibody (orb388492)

Description

Mouse monoclonal antibody to MYH11


Species/Host

Mouse

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications

Immunohistochemical
staining of human
Le...

Immunogen


Recombinant human MYH11 protein (MYH11/923);
Human uterus extract (SMMS-1)

Target

MYH11

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate
by Protein A/G. Prepared in 10mM PBS with 0.05%
rAlbumin & 0.05% azide. Also available WITHOUT
rAlbumin & azide at 1.0mg/ml.


Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody
without azide - store at -20 to -80°C. Antibody is
stable for 24 months. Non-hazardous. No MSDS
required.

Immunohistochemical
staining of human
Br...

Note

For research use only

Application notes

Smooth Muscle Myosin Heavy Chain (SM-MHC)
(Leiomyosarcoma & Myoepithelial Cell Marker)

Isotype

Mouse / IgG1, kappa + Mouse / IgG1, kappa

Clonality

Monoclonal

Immunohistochemical
staining of human
An...

MW

205kDa (MHC-1) and 200kDa (MHC-2)

Uniprot ID

P35749

Entrez

4629

Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-
P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.