

Product Datasheet

NCAM1 antibody (orb388505)

Description

Mouse monoclonal antibody to NCAM1

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Membrane preparation of a small cell lung carcinoma

Target

NCAM1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

CD56 / NCAM1 / NKH1 (Neuronal Cell Marker)

Isotype

Mouse / IgG1, kappa + Mouse / IgG1, kappa

Clonality

Monoclonal

MW

180, 145 and 125kDa

Uniprot ID

[P13592](#), [P13591](#)

Entrez


[4684](#), [24586](#)

Dilution Range


FACS: 0.5-1 µg/million cells, IF/ICC: 1-2 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human Co...



Immunohistochemical staining of human Lu...