


Product Datasheet

NEFH antibody (orb388514)

Description	Mouse monoclonal antibody to NEFH
Species/Host	Mouse
Reactivity	Gallus, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-Fr, IHC-P, WB
Immunogen	Human recombinant neurofilament protein
Target	NEFH
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	Neurofilament (NF-H) (Neuronal Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	200kDa
Uniprot ID	P12036
Entrez	4744
Dilution Range	FACS: 0.5-1 µg/million cells, IHC-P: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ce...